



Application Note 001

OpenGo Reporting Parameters

Contact

Moticon ReGo AG
Address: Machtlfinger Str. 21, 81379 Munich, Germany
EMail: support@moticon.com
Phone: +49 89 2000 301 50

Legal Note and Disclaimer

Copyright (c) 2020, Moticon ReGo AG
All rights reserved.

Any redistribution in source or binary forms, with or without modification, is not permitted.

THIS APPLICATION NOTE AND ANY REFERENCED SOFTWARE IS PROVIDED BY MOTICON REGO AG "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY, NONINFRINGEMENT, AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL MOTICON REGO AG OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

Changelog

| Version | Date | Changes |
|---------|------------|--|
| 2.0 | 01.09.2021 | Changed product name from SCIENCE to OpenGo. |
| 1.2 | 20.03.2020 | Add detailed explanation of vector parameters. |
| 1.1 | 17.10.2019 | Add spatial parameters. Add general notes on sample rates. |
| 1.0 | 14.08.2019 | Initial version |

Contents

| | | |
|----------|-----------------------|----------|
| 1 | General Notes | 4 |
| 2 | Gait Report | 4 |
| 3 | Balance Report | 7 |

1 General Notes

Coordinates Various report parameters are given in x/y coordinates, where

- x is the posterior-anterior direction, and
- y is the lateral-medial direction.

In order to identify the origin of the coordinate system, the insole outline coordinates are given in the sheet “Insole Geometry” in the gait report spreadsheet export.

Sample Rate The reporting parameters can be computed for any sample rate, however the accuracy of most parameters deteriorates with low sample rates. Note that the spatial parameters are only computed if the sample rate is at least 50 Hz (i.e. for 50 Hz and 100 Hz).

Data Channels Certain parameters require specific data channels, most notably:

- Parameters related to the pressure distribution require individual pressure values.
- All spatial parameters require angular rate and acceleration (“full” channel setup).

2 Gait Report

Table 1 — General Parameters

| Parameter | Unit | Comment |
|-------------------------------------|------|---|
| Number of steps | 1 | Sum of number of left and right steps |
| Number of steps used for statistics | 1 | If set in the report settings, excludes first/last steps of walking sequences in order to exclude acceleration/deceleration |

Table 2 — Gait Line Parameters

| Parameter | Unit | Comment |
|--|------|---|
| Mean length/width of gait line | mm | Length/width (x direction/y direction) of the average gait line |
| Mean startpoint x/y of gait line | mm | Point where the average gait line starts (typically in heel area) |
| Mean endpoint x/y of gait line | mm | Point where the average gait line ends (typically in forefoot area) |
| Standard deviation x/y of gait line start-point | mm | Indicates how much the gait line start points spread around the mean start point |
| Standard deviation x/y of gait line end-point | mm | Indicates how much the gait line end points spread around the mean end point |
| Mean gait line (anterior-posterior) ¹ | mm | x coordinates of the full representation of the average gait line as discretized points |
| Mean gait line (medial-lateral) ¹ | mm | y coordinates of the full representation of the average gait line as discretized points |
| Stddev of average gait line | mm | Standard deviation of the gait line points in y direction, given for each x coordinate of the average gait line |

¹ The *mean gait line* is determined by averaging the gait lines of individual steps, not by computing the gait line of the mean pressure distribution sequence. This way, the mean gait line can be computed in a consistent way based

on COP values, even for those sensor setups which do not include individual pressure sensors, but only the COP. The averaging procedure divides the sensor insole in 200 sections along the x axis (anterior-posterior), and determines the average COP value on the y axis (medial-lateral) for each section, separately. The length of the mean gait line is determined by the most posterior beginning and the most anterior ending of all gait lines. Numerically, the mean gait line is given as vector, and the number of elements in the vector is equal to $l_{GL}/(l_{SI}/200) + 1$, where l_{SI} is the length of the sensor insole (in mm), and l_{GL} is the length of the mean gait line (in mm). The spacing of consecutive values of the mean gait line along the x axis (anterior-posterior) is equal to $l_{SI}/200$. Consequently, the vector length is usually different for left and right. The same applies to the *standard deviation* of the mean gait line in medial-lateral direction, which is computed for each of the $l_{GL}/(l_{SI}/200) + 1$ sections, separately.

Table 3 — Pressure Distribution Parameters

| Parameter | Unit | Comment |
|-----------------------------|-------------------|---|
| Mean | N/cm ² | Mean pressure ² averaged over the entire measurement |
| Maximum | N/cm ² | Maximum pressure ² occurring during the entire measurement |
| Mean during stance phase | N/cm ² | Mean pressure ² averaged over all stance phases |
| Maximum during stance phase | N/cm ² | Maximum pressure ² occurring during all stance phases |
| During initial contact | N/cm ² | Pressure ² occurring during the initial contact, averaged over all steps |
| During mid stance | N/cm ² | Pressure ² occurring during the mid stance, averaged over all steps |
| During terminal stance | N/cm ² | Pressure ² occurring during the terminal stance, averaged over all steps |

² For each of the 16 pressure sensors.

Table 4 — Ground Reaction Force Parameters

| Parameter | Unit | Comment |
|--|------|---|
| Mean total force during stance phase | N | Mean value of the mean total force curve |
| Maximum total force during stance phase | N | Maximum of the the mean total force curve |
| Mean of all maxima of total force during all stance phases | N | This is not using the mean total force curve, but averages the maximum total force of all steps directly |
| Mean total force curve | N | Mean total force curve ³ as discretized data points, where the number of points depends on the average ground contact duration and the sample rate |
| Stddev of total force curve | N | Standard deviation of the total force curve data points over all steps, using the same duration scaling as for the mean total force curve |

³ The *mean total force curve* is the result of averaging the total force from initial contact to toe off over all detected steps. For averaging, the steps are scaled to a common duration, which is the mean stance duration t_{St} . Numerically, the mean total force curve is given as vector, and the number of elements in the vector is equal to $t_{St} \cdot f_S + 2$ (including the leading and terminating zero value), where t_{St} is the mean stance duration, and f_S is the sample rate (e.g. 100 Hz). The temporal spacing of consecutive values of the vector is equal to f_S^{-1} . The same applies to the *standard deviation* of the mean total force curve.

Table 5 — Temporal Parameters

| Parameter | Unit | Comment |
|--|-------|--|
| Mean gait cycle time | s | Time from the initial contact to the next initial contact of the same foot, averaged over all steps |
| Mean gait cadence | /min. | Mean number of strides per minute, computed using the mean gait cycle time |
| Mean double support time | s | Time during which both feet are on the ground, averaged over all steps |
| Mean fraction of double support | % | Mean double support time, relative to the mean gait cycle time |
| Mean double support time (left/right) ⁴ | s | Side-specific mean double support time, with the left mean double support time considering the time starting from the right foot's initial contact, until the left foot's toe off (and vice versa) |
| Mean step duration (left/right) ⁴ | s | Average time between consecutive initial contacts of left and right foot (and vice versa) |
| Mean stance duration | s | Time from initial contact to toe off (of the respective foot side), averaged over all steps |
| Stddev of stance duration | s | Standard deviation of the stance duration |
| Mean swing duration | s | Time from toe off to initial contact (of the respective foot side), averaged over all steps |
| Stddev of swing duration | s | Standard deviation of the swing duration |
| Mean fraction of stance phase | % | Mean stance duration, relative to the mean gait cycle time |
| Stddev of fraction of stance phase | % | Standard deviation of the mean fraction of stance phase |
| Mean fraction of swing phase | % | Mean swing duration, relative to the mean gait cycle time |
| Stddev of fraction of swing phase | % | Standard deviation of the mean fraction of swing phase |

⁴Requires high-precision synchrony between left and right data, and is therefore only applicable to short measurements.

Table 6 — Weight Bearing Parameters

| Parameter | Unit | Comment |
|--|------|--|
| Weight histogram values ⁵ | % | Histogram values representing the relative frequency of different total force ranges (bins), over the entire measurement |
| Weight histogram bins | kg | Upper border of the bins used for the weight histogram |
| Weight histogram values during stance phase ⁵ | % | Histogram values representing the relative frequency of different total force ranges (bins), limited to stance phases |
| Weight histogram bins during stance | kg | Upper border of the bins used for the weight histogram during stance phase |

⁵ The *weight histogram* uses a fixed bin width of 25 Newton to count the relative frequency of respective total force values. The number of bins is defined as 80, and loads exceeding $80 \cdot 25 \text{ N} = 200 \text{ N}$ are counted in the highest (80th) bin. Since the weight histogram is typically looked at with a notion of body weight and kilograms in mind, the bin boundaries are given in kg using a gravity rounded to 10 m/s^2 .

Table 7 — Acceleration Parameters

| Parameter | Unit | Comment |
|--|------|---|
| Mean acceleration (x/y/z) over gait cycle ⁶ | g | Result of averaging the acceleration curve over the gait cycle of all detected steps, where the gait cycles are scaled to a common duration, which is the average gait cycle time |
| Stddev of acceleration (x/y/z) over gait cycle | g | Standard deviation of the mean acceleration curve data points over all steps, using the same duration scaling as for the mean acceleration |

⁶ The *mean acceleration* over the gait cycle is the result of averaging the acceleration over the gait cycle of all detected steps. For averaging, the gait cycles are scaled to a common duration, which is the mean gait cycle time t_{GC} . Numerically, the mean acceleration curve is given as vector, and the number of elements in the vector is equal to $t_{GC} \cdot f_S$, where t_{GC} is the mean gait cycle time, and f_S is the sample rate (e.g. 100 Hz). The temporal spacing of consecutive values of the vector is equal to f_S^{-1} . The same applies to the *standard deviation* of the mean acceleration curve.

Table 8 — Spatial Parameters

| Parameter | Unit | Comment |
|--------------------|------|--|
| Mean stride length | m | The stride length is defined as the displacement of the same foot in walking direction. This parameter is the mean over all detected steps. The mean is determined from left and right data, separately, and finally averaged over left and right. |
| Walking distance | m | The walking distance traveled over the entire measurement. |
| Mean walking speed | m/s | The mean stride length divided by the mean gait cycle time. By not calculating this parameter from the walking distance and the measurement duration, non-walking periods will not affect (reduce) the mean walking speed. |

3 Balance Report

Table 9 — Center of Pressure (COP) Parameters

| Parameter | Unit | Comment |
|-----------------------------------|------|---|
| Mean COP (AP/ML) | mm | Mean center of pressure (COP) in x direction/y direction ⁷ |
| Standard deviation of COP (AP/ML) | mm | Standard deviation of the COP in x direction/y direction ⁷ |
| Bounding box of COP (AP/ML) | mm | Length/width (x direction/y direction) of the box just containing all COP points, i.e. the range of COP points ⁷ |
| Mean COP velocity | mm/s | Mean velocity of the COP travelling across the sensor insole surface ⁷ |
| COP trace length | m | Overall travel of the COP across the sensor insole surface ⁷ |

⁷ The balance report automatically detects single-leg-stance phases, and limits the computation of the parameter to these phases. If such phases are found for just one leg, the other leg is not considered. If single-leg-stance phases are not found at all, the parameter is computed over the entire measurement.