

OpenGo
by MOTICON



Sensor Insole Specification

PRELIMINARY

Effective September 1st, 2021

Content

- 03 Overview
- 04 Pressure Sensors
- 05 Pressure Sensor Geometry
- 06 Pressure Sensor Calibration
- 07 Pressure Sensor Zeroing
- 08 Acceleration Sensors
- 09 Angular Rate Sensors
- 10 Data Acquisition
- 12 Operation Modes
- 13 System Functions
- 14 Connectivity
- 16 Sizing
- 18 Power Supply
- 19 Imprint

Overview

Model
Insole3



Moticon's *OpenGo Sensor Insole* is a versatile tool for sensing human foot dynamics.


The sensor insole is fully integrated and highly flexible. It incorporates all technological components to execute standalone measurements. No external devices are required for data acquisition.

The following table depicts the basic specifications. For more detailed information, refer to the corresponding chapters of this document.

Sensors	16 pressure 3 acceleration 3 angular rate per side
Wireless	Bluetooth Low Energy
Power Supply	PD2032 coin cell battery, rechargeable
Sizes	9 double sizes EU 32/33 - 48/49 US 1/2 - 12½/13½
Data Storage	On-board memory & live transmission to smartphone

Pressure Sensors

Quantity	Type	Axis	Range	Resolution	Hysteresis	Calibration
16 per side	capacitive plantar	Z	0-50 N/cm ²	0.25 N/cm ²	≤ 1 %	in-production lifetime



Notes

The sensor insole outline and the sensor positions are given in a common coordinate system. The orientation of the coordinate systems is different for left and right sides.

Legend

- ① Sensor numbering
- X/Y/Z Directions of coordinate system
- CL/R Center of coordinate system

Pressure Sensor Geometry

	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8	Size 9	
TAI	12,129	13,240	14,462	15,787	17,237	18,824	20,572	22,456	24,524	mm ²
TAS	7,578	8,330	9,179	10,109	11,129	12,241	13,696	15,042	16,525	mm ²
SC ¹	62.5	62.9	63.5	64.0	64.6	65.0	66.6	67.0	67.4	%

Notes

¹ Sensor coverage varies over size due to keepout areas for wiring, electronics, antenna and borders which do not scale with size.

Legend

TAI: Total area of sensor insole.

TAS: Total area covered by pressure sensors.

SC: Sensor coverage. Area covered by pressure sensors in relation to total area of sensor insole.

Pressure Sensor Calibration


Type	Device	Procedure
in-production testing (IPT)	fluid cushion	each sensor insole is calibrated individually
user calibration	n/a	user to carry out motion sequence to calibrate sensor insole individually for highest accuracy ¹

Notes

¹The calibration routine is carried out using the *OpenGo App*. Visit moticon.com/doc-opengo-calibration for more information.

Example

The figure shows a typical stack of calibration curves as generated in the in-production testing, representing the pressure sensor response function when loaded. Each calibration curve relates to one discrete pressure sensor.



Pressure Sensor Zeroing

Type	Description
automatic zeroing mode ¹	offsets and drifts in the pressure sensor readings are compensated automatically by software routines.
manual zeroing ²	manual zeroing of pressure sensor readings.

Notes

¹ The automatic zeroing mode is always active. It is based on algorithms continuously checking sensor zero levels and compensates for sensor offsets and drifts which may occur due to lacing shoes and due to temperature changes.


² In addition, a manual zeroing routine can be carried out by means of a software function in the *OpenGo App*. The routine does not disable the automatic zeroing, but overwrites current zero values. The routine should be applied in cases where no unloading cycles occur or are not feasible before measurement start.

Example

The following example shows the effect of the automatic zeroing mode.

The graph shows the total force readings of the left (blue) and right (brown) leg of a normal gait.

After 2-5 steps of walking, initial offsets are compensated to a level between 0N and 5N which is caused by sensor noise only.



Acceleration Sensors

Quantity	Type	Axis	Range	Resolution	Hysteresis	Calibration
1 per side	MEMS LSM6DSL	X/Y/Z	± 16 g	0.488 mg/LSB 16 bit	n.a.	in-factory lifetime

Notes

The position of the IMU sensor is at the origin of the coordinate system.

The orientation of the coordinate system is different for left and right sides.

Any coordinates to the sensor centers and outline coordinates of the data output are given in relation to the coordinate system.


Refer to the official LSM6DSL ST Microelectronics datasheet for further information.

Legend

① Sensor numbering

X/Y/Z Directions of coordinate system and acceleration

C_{L/R} Center of coordinate system



Angular Rate Sensors

Quantity	Type	Axis	Range	Resolution	Hysteresis	Calibration
1 per side	MEMS LSM6DSL	$\Omega_X/\Omega_Y/\Omega_Z$	± 2000 dps	70 mdsp/LSB 16 bit	n.a.	in-factory lifetime


Notes

The position of the IMU sensor is at the origin of the coordinate system.
The orientation of the coordinate system is different for left and right sides.
Any coordinates to the sensor centers and outline coordinates of the data output are given in relation to the coordinate system.

Refer to the official LSM6DSL ST Microelectronics datasheet for further information.

Legend

- ① Sensor numbering
- X/Y/Z Directions of coordinate system and acceleration
- $\Omega_{X/Y/Z}$ Angular rates rotation directions
- $C_{L/R}$ Center of coordinate system



Data Acquisition

Sample Rates ¹	Memory	Timer ²	Primary Data Channels ¹	Computed Data Channels ¹
10 Hz	32 MB ³	drift ≤ 1 %	TS ²	COP [X,Y]
25 Hz	FLASH		PR [1...16]	TF [Z]
50 Hz	non-		ACC [X,Y,Z]	
100 Hz	extendable		ANG [X,Y,Z]	

Notes

¹ Sample rates and data channel combinations are software-selectable in the *OpenGo App*.

Each data channel can be activated/deactivated individually. The time stamp channel is always active.

² Left and right sensor insoles are time synchronized. Time drift refers to total drift of timing data over time from start of measurement.

Date information (yyyy/mm/dd) from control device, i.e. mobile device/ smartphone, is recorded upon start of measurement.

For the time stamp, the internal clock creates timing information at the resolution of the sample rate (hh:mm:ss.xyz).

³ From production date 08-2020. Before that date: 16 MB.

Legend

TS: Time stamp in s (resolution in ms).

PR: Pressure channels 1 to 16 in N/cm².

ACC: Acceleration channels in X, Y and Z axis direction in g.

ANG: Angular rates about the X, Y and Z axes in dps.

COP: Center of pressure in X and Y axis direction in the range [-0.5...+0.5] related to corresponding sensor insole length or width.

TF: Total force (= ground reaction force) in Z axis direction in N.

Data Acquisition

Sensor Setup ¹	TR @ 10 Hz	TR @ 25 Hz	TR @ 50 Hz	TR @ 100 Hz	
Full Setup	11:20	04:32	02:16	01:08	hh:mm
Standard Setup	13:58	05:35	02:47	01:24	hh:mm
Monitoring Setup	32:15	12:54	06:27	03:13	hh:mm

Sensor Setup ¹	TS	PR	ACC	ANG	COP	TF
Full Setup	●	●	●	●	●	●
Standard Setup	●	●	●		●	●
Monitoring Setup	●		● ²		●	●

Notes

Recording times refer to continuous data recording in Recording mode. Effective recording times can be higher when using SmartRecording mode for activity triggered data recording. For example, A measurement acquired with standard sensor setup in SmartRecording mode @ 50 Hz can take a up to 5 times of the duration of continuous recording if the wearer of the sensor insoles is only active 20 % of the time.

¹ The depicted sensor setups are setups which are pre-defined in the *OpenGo App*.

Other Sample rates and data channel combinations are software-selectable in the *OpenGo App*.

Each data channel can be activated/deactivated individually. The time stamp channel (TM) is always active.

Indicated recording times apply for sensor insoles with 32 MB memory produced from 08/2020!

² Only Z axis.

Legend

Tr: Maximum recording time for continuous on-board data recording.

Operation Modes

Standard Modes ¹	Description	Use Case
Preview	Transmits fixed basic sensor data wirelessly	Live showcasing and activity checks in the <i>OpenGo App</i>
Live	Transmits sensor data wirelessly to endpoint (cloud/desktop computer)	Direct data transmission for storage on endpoint
Record	On-board sensor data recording	Standalone data acquisition without smartphone
Transmit	Wireless sensor data transmission from on-board memory to endpoint PC	Transmitting previously recorded sensor data to endpoint for storage on endpoint

Intelligence Modes ¹	Description	Use Case
Smart Recording ²	Activity triggered recording mode. Recording is only on when users are active.	For longterm gait monitoring
Smart Sleep	Automatic shut down into power safe mode when sensor insole is not in use. No hardware switches are used.	No user interaction for turning on/off required

Notes

¹ The operation modes can be selected in the *OpenGo App*.

² Smart Recording is an activity triggered recording mode where internal firmware algorithms control the intermediate pausing and resuming in action. The mobile device (phone/tablet) is only required for starting and stopping, not for the recording itself. For further details on the operation modes refer to the product pages on moticon.com/opengo/app.

System Functions

Function	Description
Firmware Update	Firmware of the sensor insoles can be updated wirelessly using the <i>OpenGo App</i> .
File Management	The firmware runs a file management system on the on-board memory of the sensor insoles. Measurements can be selected, deleted and transmitted individually using the <i>OpenGo App</i> .
System Status	The firmware checks the system status and transmits status information. Status information includes: <ol style="list-style-type: none">1. Restarts2. Active time3. Miscellaneous: firmware version, serial number, size, side
Battery Status	The firmware detects the battery charge status of the inserted coin cell batteries. The charge status is transmitted to the <i>OpenGo App</i> .

Connectivity

Wireless Standard	Profile	Wireless Range	Bandwidth
BLE 5.0 ¹	proprietary ²	≥ 5.0 m ³	54 kB/s ⁴

Data Volume	Download Time ⁵	Recording Equivalent ⁵
28.80 MB ⁶	00:08:54 hh:mm:ss	68 min @ 100 Hz Full Setup
4.30 MB	00:01:20 hh:mm:ss	10 min @ 100 Hz Full Setup
0.43 MB	00:00:08 hh:mm:ss	1 min @ 100 Hz Full Setup

Notes

¹ Version refers to integrated wireless chip. BLE connection to smartphones can be established from version 5.0.

² No standard BLE profile is used.

³ Depends on environmental conditions and other active wireless devices.

⁴ Maximum data transmission bandwidth. Depends on smartphone performance.

⁵ Given times and setups are estimates and shall illustrate expectable download times of previously on-board recorded sensor data to a smartphone using the *OpenGo App*.

⁶ Includes 10 % overhead for auxiliary data.

Connectivity



Notes


For instructions of use refer to the tutorials on moticon.com/documentation.

Legend

--- Wireless data transmission.

Sizes

	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8	Size 9
EU	32/33	34/35	36/37	38/39	40/41	42/43	44/45	46/47	48/49
US	1/2	2½/3	3½/4½	5/6	6½/7½	8/9	9½/10½	11/12	12½/13½




Notes

Refer to the *OpenGo Sensor Insole* fitting guide for more detailed information moticon.com/opengo/sensor-insole-sizes.

Sizes

	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8	Size 9	
W	80.2	83.4	86.7	90.2	93.8	97.5	101.4	105.5	109.7	mm
L	214.9	225.6	236.8	248.6	261.1	274.2	288.0	302.4	317.5	mm
T ₁	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	mm
T ₂	10.5	10.5	10.5	10.5	10.5	10.5	10.5	10.5	10.5	mm




Notes

Dimensions may vary ± 1.5 mm due to production tolerances.

Power Supply

Supply Voltage	Power Consumption	Battery ¹ Type	Battery Form Factor	Battery Technology
3.7 v ± 0.4 v	n/a mAh	1 coin cell rechargeable	2032	Lithium Ion



Caution:

¹ Recommended coin cell battery type is PD2032 (Route Jade, South Korea), available from Moticon.
DO NOT USE primary coin cell batteries type CR2032 or mangan based rechargeable coin cell batteries type ML2032.
These batteries have a nominal voltage of 3.0 V and continuous voltage drop in use.
Moticon SCIENCE Sensor Insoles require a minimum operation voltage of 3.2 V. A lower operating voltage causes malfunction.

Notes

Refer to the tutorials for coin cell battery handling and charging recommendations on moticon.com/opengo/documentation.

Imprint

Legal Notice

This document is the property of Moticon ReGo AG.
Any reproduction of this document in part or in whole is strictly prohibited.
For educational purposes only.
Errors and omissions excepted and subject to change without notice.

Further Information

For more detailed product information, please visit
moticon.com/opengo/sensor-insole
moticon.com/opengo/resources
moticon.com/opengo/documentation

Release Notes

Version 03.01
Date September 1st 2021

Vendor Information

Moticon ReGo AG
Machtlfinger Str. 21
81379 Munich
Germany

Phone +49 89 2000 301 60
EMail sales@moticon.com
Web moticon.com



Moticon
Sole
IP X1
Model: Insole 3
FCC ID: 687316-ansole3
OpenGo - RIGHT
6 - 42/43 - 8/9



moticon

Email sales@moticon.com
Web moticon.com/opengo