

**OpenGo**  
by MOTICON



# Price List

Germany & EU

Effective September 1st, 2021

# Content

- 03 Sensor Insoles  
Spare Parts
- 04 Power Supply  
Product Bags
- 05 Desktop Software  
Mobile App
- 06 Services
- 07 Imprint

# Sensor Insoles

Product	Description	PU	Price (Euro)
OpenGo Sensor Insoles	Connectivity Bluetooth LE 4.2 or higher Sample Rates 10, 25, 50, 100 Hz Sensors 16 pressure sensors 6-axis IMU Memory 32MB internal Data Acquisition Recording, SmartRecording, Live, SmartSleep Data Output Pressure, total force, COP, acceleration, angular rates, time Power Supply Li-ion 2032 coin cell battery, 3.7V Sizes 9 double sizes	1 pair	1495.00

# Spare Parts

Product	Description	PU	Price (Euro)
OpenGo Battery Lid	Replacement, fits OpenGo Sensor Insoles	1 pair	20.00

# Power Supply

Product	Description	PU	Price (Euro)
Coin Cell Battery	PD2032, Li-ion, rechargeable Voltage 3.7V, capacity 75 mAh	1 pair	20.00
Moticon Coin Cell Charger	2 charging slots for PD2032 coin cells Charge voltage 4.2V, Supply voltage 5.0V Micro-USB Dimensions 46x34x23 mm LED charge status indicator	1 piece	95.00

# Product Bags

Product	Description	PU	Price (Euro)
OpenGo Product Bag	Compartments for 6 pairs of sensor insoles 2 flat compartments, 1 large side compartment Outer material Cordura, carry strap, zippers Dimensions 430x300x60 mm	1 piece	145.00
OpenGo Coin Cell Bag	Compartments for 20 PD2032 coin cells Outer material Cordura, zippers Dimensions 150 x 90 x 15 mm	1 piece	25.00
OpenGo Sensor Insole Sack	Fits 1 pair of sensor insoles Outer material mesh textile Dimensions 310 x 120 mm	1 piece	15.00

# Software

Product	Description	PU	Price (Euro)
OpenGo Desktop Software <sup>1</sup>	Basic data management & analysis	1 license	995.00
OpenGo Analyze Module <sup>2</sup>	Advanced data labeling, analysis & processing	1 license	995.00
OpenGo Video Module <sup>2</sup>	Automated video synchronization Multi video channel management	1 license	995.00
OpenGo Gait Report <sup>2</sup>	Automated gait report generation Custom report setup & export functions	1 license	495.00
OpenGo Jump Report <sup>2</sup>	Automated jump report generation Custom report setup & export functions	1 license	495.00
OpenGo Balance Report <sup>2</sup>	Automated balance report generation Custom report setup & export functions	1 license	495.00

<sup>1</sup> Minimum system requirement: Windows 10, Intel i5 or equivalent, 4 GB RAM, 1 GB HDD, DirectX9

<sup>2</sup> Requires OpenGo Desktop Software license

# App

Product	Description	PU	Price (Euro)
OpenGo Mobile App <sup>3</sup>	Sensor insole control and settings Measurement organization, data transfer to OpenGo Software	1 license	0.00

<sup>3</sup> For Android mobile devices only, available in Google PlayStore for multiple accounts, registration required  
System requirements: Android 9 (Pie), BLE 5.0

# Services

<b>Product</b>	<b>Description</b>	<b>Unit</b>	<b>Price (Euro)</b>
Support Service	12 months customer support Mandatory for 1st order	1 service	295.00
Software Update Service	12 months updates for OpenGo App & Software Mandatory for 1st order	1 service	295.00
Moticon Sensor Insole Replacement Service	On-demand service linked to purchase of each pair of Moticon Sensor Insoles Valid 12 months after purchase, unlimited exchanges	1 service per pair	1195.00

# Imprint

## Price List Germany & EU

### Legal Notice

This document is the property of Moticon ReGo AG. Any reproduction of this document in part or in whole is strictly prohibited. For educational purposes only. Errors and omissions excepted and subject to change without notice.

### Notes on Prices

This price list informs you about delivery ex works (EXW).

Prices are net plus VAT.

For international orders, the VAT must be paid locally by the customer.

For non-EU orders, customs duty fees may apply according to national import tax rates.

The Moticon general terms and conditions (GTC) under [www.moticon.com/terms-and-conditions](http://www.moticon.com/terms-and-conditions) apply to this price list.

Moticon offers international carriage, insurance and custom services, but cannot be held liable for any loss that may occur during the carriage.

### Further Information

For more detailed product information, please visit [www.moticon.com/opengo](http://www.moticon.com/opengo)

### Release Notes

Version 03.01

### Vendor Information

Moticon ReGo AG  
Machtlfinger Str. 21  
81379 Munich  
Germany

Phone **+49 89 2000 301 60**  
EMail **[sales@moticon.com](mailto:sales@moticon.com)**  
Web **[www.moticon.com](http://www.moticon.com)**



Moticon  
Germany  
2017/11/14

CE

SN 6589

Model: Innote3 IP X1

FCC ID: 687316-Innote3

OpenGo - RIGHT

6 - 42/43 - 8/9



MOTICON

Email [sales@moticon.com](mailto:sales@moticon.com)  
Web [moticon.com/opengo](http://moticon.com/opengo)