



OpenGo
byMOTICON

Price List

United States

Effective September 1st, 2022

Content

03 Sensor Insoles
Spare Parts

04 Power Supply
Product Bags

05 Desktop Software
Mobile App

06 Services

Sensor Insoles

Product	Description	PU	Price (Dollar)	
OpenGo Sensor Insoles	Connectivity	Bluetooth LE 4.2 or higher	1 pair	1,995.00
	Sample Rates	10, 25, 50, 100 Hz		
	Sensors	16 pressure sensors 6-axis IMU		
	Memory	32MB internal		
	Data Acquisition	Recording, SmartRecording, Live, SmartSleep		
	Data Output	Pressure, total force, COP, acceleration, angular rates, time		
	Power Supply	Li-ion 2032 coin cell battery, 3.7V		
	Sizes	9 different double sizes available		

Spare Parts

Product	Description	PU	Price (Dollar)
OpenGo Battery Lid	Replacement, fits OpenGo Sensor Insoles	1 pair	20.00

Power Supply

Product	Description	PU	Price (Dollar)
Coin Cell Battery	PD2032, Li-ion, rechargeable Voltage 3.7 V, capacity 75 mAh	1 pair	25.00
Moticon Coin Cell Charger	2 charging slots for PD2032 coin cells Supply voltage 5.0 V Micro-USB Charge voltage 4.2 V Dimensions 46 x 34 x 23 mm LED charge status indicator	1 piece	115.00

Product Bags

Product	Description	PU	Price (Dollar)
OpenGo Product Bag	Compartments for 6 pairs of sensor insoles 2 flat compartments, 1 large side compartment Outer material Cordura, carry strap, zippers Dimensions 430 x 300 x 60 mm	1 piece	195.00
OpenGo Coin Cell Bag	Compartments for 20 PD2032 coin cells Outer material Cordura, zippers Dimensions 150 x 90 x 15 mm	1 piece	25.00
OpenGo Sensor Insole Sack	Fits 1 pair of sensor insoles Outer material mesh textile Dimensions 310 x 120 mm	1 piece	15.00

Software

Product	Description	PU	Price (Dollar)
OpenGo Desktop Software ¹	Basic data management & analysis	1 license	1,195.00
OpenGo Analyze Module ²	Advanced data labeling, analysis & processing	1 license	1,195.00
OpenGo Video Module ²	Automated video synchronization Multi video channel management	1 license	1,195.00
OpenGo Gait Report ²	Automated gait report generation Custom report setup & export functions	1 license	595.00
OpenGo Jump Report ²	Automated jump report generation Custom report setup & export functions	1 license	595.00
OpenGo Balance Report ²	Automated balance report generation Custom report setup & export functions	1 license	595.00

¹ Minimum system requirement: Windows 10, Intel i5 or equivalent, 4 GB RAM, 1 GB HDD, DirectX9

² Requires OpenGo Desktop Software license

App

Product	Description	PU	Price (Dollar)
OpenGo Mobile App ³	Sensor insole control and settings Measurement organization Data transfer to OpenGo Software	1 license	0.00

³ For Android mobile devices only, available in Google PlayStore for multiple accounts, registration required
System requirements: Android 9 (Pie), BLE 5.0

Services

Product	Description	Unit	Price (Dollar)
Support Service	12 months customer support (mandatory for initial order)	1 service	395.00
Software Update Service	12 months updates for OpenGo App & Software (mandatory for initial order)	1 service	395.00
Moticon Sensor Insole Replacement Service	On-demand service for product exchange (not required for initial order), includes shipment and return of exchanged product	1 service per pair	1595.00

Legal Notes

Errors and omissions excepted .

Prices may change without prior notice.

Prices are net prices and local VAT applies.

The Moticon General Terms and Conditions apply moticon.com/terms-and-conditions.

Product Information

moticon.com/opengo

Sensor Insole Information

moticon.com/opengo/sensor-insoles

Software Information

moticon.com/opengo/software

App Information

moticon.com/opengo/app

Accessories Information

moticon.com/opengo/accessories

Services Information

moticon.com/opengo/services



Original English Version
Keep this document confidential

Moticon Rego AG
Machtlfinger Str. 21
81379 Munich
Germany

+49 89 200 030 160
sales@moticon.com
moticon.com/opengo

REV 03.02