



Sensor Insole

Fitting Guide

Effective June 1st, 2025

Content

- 03 Use Notes
- 04 Size Overview
- 05 Sensor Insole Printouts
- 15 Legal Notes & Links

Use Notes


This document enables you to check the sizes of Moticon OpenGo Sensor Insoles by printing the outlines on paper in real size.

Printing

Make sure to **print** the insole outline **without scaling** to achieve proper 1:1 sizes. It is recommended to print on heavy paper ($\geq 120\text{g}/\text{sheet}$) if you intend to insert the template into a shoe.

Cutting


Cut with scissors or paper knife along the outer line in order to achieve the right outline shape. Make sure you cut or fold the dividing line between the two halves to achieve the right length.



Size Overview


Size 1

EU 32 - 33 US Men 1 - 2 US Women 2½ - 3½




Size 2

EU 34 - 35 US Men 2½ - 3 US Women 4 - 4½



Size 3

EU 36 - 37 US Men 3½ - 4½ US Women 5 - 6




Size 4

EU 38 - 39 US Men 5 - 6 US Women 6½ - 7½




Size 5

EU 40 - 41 US Men 6½ - 7½ US Women 8 - 9




Size 6

EU 42 - 43 US Men 8 - 9 US Women 9½ - 10½




Size 7

EU 44 - 45 US Men 9½ - 10½ US Women 11 - 12




Size 8



EU 46 - 47 US Men 11 - 12 US Women 12½ - 12½




Size 9

EU 48 - 49 US Men 12½ - 13½ US Women 14 - 15







Size 1

EU 32 - 33

US Men 1 - 2

US Women 2½ - 3½

Width


80.2 mm

3.16 in

Length

214.9 mm

8.46 in



Size 2

EU 34 - 35

US Men 2½ - 3

US Women 4 - 4½

Width



83.4 mm


3.28 in

Length

225.6 mm

8.88 in





Size 3

EU 36 - 37

US Men 3½ - 4½

US Women 5 - 6

Width


86.7 mm

3.41 in

Length

236.8 mm

9.32 in



Size 4

EU 38 - 39

US Men 5 - 6

US Women 6½ - 7½

Width



90.2 mm


3.55 in

Length

248.6 mm

9.79 in





Size 5

EU 40 - 41

US Men 6½ - 7½

US Women 8 - 9

Width


93.8 mm

3.69 in

Length

261.1 mm

10.28 in



Size 6

EU 42 - 43

US Men 8 - 9

US Women 9½ - 10½

Width


97.5 mm

3.84 in


Length

274.2 mm


10.80 in



OpenGo
by MOTICOR



OpenGo
by MOTICOR



Size 7

EU 44 - 45

US Men 9½ - 10½

US Women 11 - 12

Width


101.4 mm

3.99 in

Length

288.0 mm

11.34 in



Size 8

EU 46 - 47

US Men 11 - 12

US Women 12½ - 12½

Width


105.5 mm


4.15 in

Length

302.4 mm

11.91 in





Size 7

EU 48 - 49

US Men 12½ - 13½

US Women 14 - 15

Width

109.7 mm

4.32 in

Length

317.5 mm

12.50 in

Legal Notes

Errors and omissions excepted.

Product Information

moticon.com/opengo

Sensor Insole Information

moticon.com/opengo/sensor-insoles

Software Information

moticon.com/opengo/software

App Information

moticon.com/opengo/app

Accessories Information

moticon.com/opengo/accessories

Services Information

moticon.com/opengo/services



Original English Version
June 2025

Moticon ReGo AG
Machtlfinger Str. 21
81379 Munich
Germany

+49 89 200 030 160
sales@moticon.com

REV 03.02