

**OpenGo**  
by MOTICON

# Sensor Insole Fitting Guide

Effective September 1st, 2021



# Content

03 Size Overview

04 Size 1 Details

06 Size 2 Details

08 Size 3 Details

10 Size 4 Details

12 Size 5 Details

14 Size 6 Details

16 Size 7 Details

18 Size 8 Details


20 Size 9 Details

22 Use Notes

23 Imprint


# Size Overview


	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8	Size 9
EU	32/33	34/35	36/37	38/39	40/41	42/43	44/45	46/47	48/49
US	1/2	2½/3	3½/4½	5/6	6½/7½	8/9	9½/10½	11/12	12½/13½



# Size 1


	Size		Width	Length
EU	32/33	mm	80.2	214.9
US	1/2	in	3.16	8.46






# Size 2


	Size		Width	Length
EU	34/35	mm	83.4	225.6
US	2½/3	in	3.28	8.88






# Size 3

	Size		Width	Length
EU	36/37	mm	86.7	236.8
US	3½/4½	in	3.41	9.32









# Size 4


	Size		Width	Length
EU	38/39	mm	90.2	248.6
US	5/6	in	3.55	9.79






# Size 5


	Size		Width	Length
EU	40/41	mm	93.8	261.1
US	6½/7½	in	3.69	10.28






# Size 6


	Size		Width	Length
EU	42/43	mm	97.50	274.2
US	8/9	in	3.84	10.80






# Size 7

	Size		Width	Length
EU	44/45	mm	101.4	288.0
US	9½/10½	in	3.99	11.34









# Size 8


	Size		Width	Length
EU	46/47	mm	105.5	302.4
US	11/12	in	4.15	11.91






# Size 9

	Size		Width	Length
EU	48/49	mm	109.7	317.5
US	12½/13½	in	4.32	12.50





# Use Notes

## General


This document enables you to check the sizes of Moticon SCIENCE Sensor Insoles by printing the outlines on paper in real size.

## Printing

Make sure to **print** the insole outlines **without scaling** to achieve proper 1:1 sizes. It is recommended to print on heavy paper ( $\geq 120\text{g/sheet}$ ) if you intend to insert the template into a shoe.

## Cutting

Cut with scissors or paper knife along the outer line in order to achieve the right outline shape. Make sure you cut or fold the dividing line between the two halves to achieve the right length.



# Imprint

## Legal Notice

This document is the property of Moticon ReGo AG.  
Any reproduction of this document in part or in whole is strictly prohibited.  
For educational purposes only.  
Errors and omissions excepted and subject to change without notice.

## Further Information

For more detailed product information, please visit  
[moticon.com/opengo](https://moticon.com/opengo)  
[moticon.com/opengo/documentation](https://moticon.com/opengo/documentation)  
[moticon.com/opengo/resources](https://moticon.com/opengo/resources)

## Release Notes

Version 03.01  
Date September 1st, 2021

## Vendor Information

Moticon ReGo AG  
Machtlfinger Str. 21  
81379 Munich  
Germany

Phone +49 89 2000 301 60  
EMail [sales@moticon.com](mailto:sales@moticon.com)  
Web [moticon.com](https://moticon.com)



Email [sales@moticon.com](mailto:sales@moticon.com)  
Web [moticon.com](http://moticon.com)