



# Sensor Insoles Fitting Guide

Effective January 4th, 2021

# Use Notes

This document enables you to check the sizes of Moticon ReGo Sensor Insoles by printing the outlines on paper in real size.


## Printing

Make sure to **print** the insole outlines **without scaling** to achieve proper 1:1 sizes.

It is recommended to print on heavy paper (120g/sheet) if you intend to insert the template into a shoe.

## Cutting

Cut with scissors or paper knife along the outer line in order to achieve the right outline shape. Make sure you cut or fold the dividing line between the two halves to achieve the right length.




# ReGo Size Overview

---

## SIZE 1

EU 32 – 33 US Men 1 – 2 US Women 2½ – 3½



---

## SIZE 2

EU 34 – 35 US Men 2½ – 3 US Women 4 – 4½



---

## SIZE 3

EU 36 – 37 US Men 3½ – 4½ US Women 5 – 6



---

## SIZE 4


EU 38 – 39 US Men 5 – 6 US Women 6½ – 7½



---

## SIZE 5

EU 40 – 41 US Men 6½ – 7½ US Women 8 – 9



---

## SIZE 6


EU 42 – 43 US Men 8 – 9 US Women 9½ – 10½



---

## SIZE 7


EU 44 – 45 US Men 9½ – 10½ US Women 11 – 12



---

## SIZE 8

EU 46 – 47 US Men 11 – 12 US Women 12½ – 13½






---

## SIZE 9

EU 48 – 49 US Men 12½ – 13½ US Women 14 – 15







## ReGo 1

EU 32 – 33

US Men 1 – 2


US Women 2½ – 3½

**WIDTH**

80.2 mm  
3.16 in

**LENGTH**

214.9 mm  
8.46 in



## ReGo 2

EU 34 – 35

US Men 2½ – 3



US Women 4 – 4½


**WIDTH**

83.4 mm  
3.28 in

**LENGTH**

225.6 mm  
8.88 in






## ReGo 3

**EU 36 – 37**

**US Men 3½ – 4½**

**US Women 5 – 6**

<b>WIDTH</b>	<b>LENGTH</b>
86.7 mm 3.41 in	236.8 mm 9.32 in





## ReGo 4

**EU 38 – 39**


**US Men 5 – 6**

**US Women 6½ – 7½**

<b>WIDTH</b>	<b>LENGTH</b>
90.2 mm 3.55 in	248.6 mm 9.79 in







## ReGo 5

**EU** 40 – 41

**US Men** 6½ – 7½

**US Women** 8 – 9

### WIDTH


93.8 mm

3.69 in

### LENGTH

261.1 mm

10.28 in



## ReGo 6

**EU** 42 – 43

**US Men** 8 – 9

**US Women** 9½ – 10½

### WIDTH



97.50 mm


3.84 in

### LENGTH

274.2 mm

10.80 in





## ReGo 7

**EU 44 – 45**

**US Men 9½ – 10½**


**US Women 11 – 12**

### WIDTH

101.4 mm  
3.99 in

### LENGTH

288.0 mm  
11.34 in



## ReGo 8

**EU 46 – 47**

**US Men 11 – 12**


**US Women 12½ – 13½**


### WIDTH

105.5 mm  
4.15 in

### LENGTH

302.4 mm  
11.91 in





## ReGo 9

**EU 48 – 49**

**US Men 12½ – 13½**

**US Women 14 – 15**

---

### **WIDTH**

109.7 mm  
4.32 in

### **LENGTH**

317.5 mm  
12.50 in



Original English Version  
Keep this document confidential  
January 2021

Moticon ReGo AG  
Machtlfinger Str. 21  
81379 Munich  
Germany

+49 89 200 030 160  
sales@moticon.de  
moticon.de/support

REV 01.02