



Sensor Insoles Fitting Guide

Effective January 4th, 2021

Use Notes

This document enables you to check the sizes of Moticon ReGo Sensor Insoles by printing the outlines on paper in real size.


Printing

Make sure to **print** the insole outlines **without scaling** to achieve proper 1:1 sizes.

It is recommended to print on heavy paper (120g/sheet) if you intend to insert the template into a shoe.

Cutting


Cut with scissors or paper knife along the outer line in order to achieve the right outline shape. Make sure you cut or fold the dividing line between the two halves to achieve the right length.



ReGo Size Overview

SIZE 1

EU 32 – 33 US Men 1 – 2 US Women 2½ – 3½



SIZE 2

EU 34 – 35 US Men 2½ – 3 US Women 4 – 4½



SIZE 3

EU 36 – 37 US Men 3½ – 4½ US Women 5 – 6



SIZE 4

EU 38 – 39 US Men 5 – 6 US Women 6½ – 7½



SIZE 5

EU 40 – 41 US Men 6½ – 7½ US Women 8 – 9



SIZE 6

EU 42 – 43 US Men 8 – 9 US Women 9½ – 10½



SIZE 7

EU 44 – 45 US Men 9½ – 10½ US Women 11 – 12



SIZE 8



EU 46 – 47 US Men 11 – 12 US Women 12½ – 13½




SIZE 9

EU 48 – 49 US Men 12½ – 13½ US Women 14 – 15







ReGo 1

EU 32 – 33

US Men 1 – 2

US Women 2½ – 3½

WIDTH


80.2 mm

3.16 in

LENGTH

214.9 mm

8.46 in



ReGo 2

EU 34 – 35

US Men 2½ – 3

US Women 4 – 4½

WIDTH



83.4 mm


3.28 in

LENGTH

225.6 mm

8.88 in






ReGo 3

EU 36 – 37

US Men 3½ – 4½

US Women 5 – 6

WIDTH	LENGTH
86.7 mm 3.41 in	236.8 mm 9.32 in





ReGo 4


EU 38 – 39

US Men 5 – 6

US Women 6½ – 7½

WIDTH	LENGTH
90.2 mm 3.55 in	248.6 mm 9.79 in






ReGo 5

EU 40 – 41

US Men 6½ – 7½

US Women 8 – 9

WIDTH	LENGTH
93.8 mm	261.1 mm
3.69 in	10.28 in





ReGo 6


EU 42 – 43

US Men 8 – 9

US Women 9½ – 10½

WIDTH	LENGTH
97.50 mm	274.2 mm
3.84 in	10.80 in





ReGo 7

EU 44 – 45

US Men 9½ – 10½


US Women 11 – 12

WIDTH

101.4 mm
3.99 in

LENGTH

288.0 mm
11.34 in



ReGo 8

EU 46 – 47

US Men 11 – 12


US Women 12½ – 13½


WIDTH

105.5 mm
4.15 in

LENGTH

302.4 mm
11.91 in





ReGo 9

EU 48 – 49

US Men 12½ – 13½

US Women 14 – 15

WIDTH

109.7 mm

4.32 in

LENGTH

317.5 mm

12.50 in



Original English Version
Keep this document confidential
January 2021

Moticon ReGo AG
Machtlfinger Str. 21
81379 Munich
Germany

+49 89 200 030 160
sales@moticon.de
moticon.de/support

REV 01.02