






# Sensor Insoles

**599 \$**

*per pair*

Sensors	16 capacitive pressure sensors, 6-axis inertial measurement unit
Sample rate	100 Hz
Data acquisition	Live stream to mobile device
Connectivity	Bluetooth Low Energy 5.0
Power supply	Li-ion PD2032 coin cell battery, 3.7 V, rechargeable, USB-C connector
Sizes	9 double sizes S1 - S9 (EU 32 - 49; US men 1 - 13.5)



# Tests

**0.99 \$**

*per test*

Detailed Test Report.  
Test execution and biofeedback free of charge.

Available tests	<a href="https://moticon.com/rego/tests">moticon.com/rego/tests</a>
Free tests	10 upon initial order
Quota	No minimum quota, no subscription

# Accessories & Services

## 1 POWER MANAGEMENT

Power Supply	1 pc	Power supply 240/110 V DC / 5 V AC, USB 2.0 type A, fits 4 dual charge cables	99.00 \$
Dual Charge Cable	1 pc	USB 2.0 type A to dual type C charge cable for 1 pair of sensor insoles	9.90 \$

## 2 BAGS & MARKERS

Product Bag	1 pc	Cordura zipper bag for 6 pairs of sensor insoles & accessories	149.00 \$
Sensor Insole Sack	1 pc	Mesh sack with drawcord lock, for 1 pair of sensor insole, fits all sizes	9.90 \$
Floor Marker Disc	1 pair	Flexible marker with rubberized back, indicates clockwise & counter clockwise, diameter 18 cm / 7.1 in	19.90 \$
Floor Marker Square	1 pc	Flexible marker with rubberized back, indicates jump area, square 40 x 40 cm / 15.7 x 15.7 in	29.90 \$
Tape Measure	1 pc	For setting out gait corridors	9.90 \$

## 3 CONSUMABLES

Test Socks	100 pcs	Disposable socks, antibacterial, one size fits all	29.00 \$
Disinfectant Spray	150 ml	Sterillium disinfectant spray	9.90 \$
Disinfectant Refill	1 L	Sterillium disinfectant	29.90 \$

## 4 SERVICES

ReGo Setup	1 setup	Initial customer setup, account activation	99.00 \$
ReGo Product Training	2 hr	Comprehensive online training and Q&A, covers mobile app functions, sensor insoles and maintenance	299.00 \$



Original English Version  
Keep this document confidential  
January 2023

Moticon ReGo AG  
Machtlfinger Str. 21  
81379 Munich  
Germany

+49 89 200 030 160  
sales@moticon.com  
moticon.com/rego

Document valid for: Europe  
REV 01.05