



**ReGo**  
by MOTICON

# Price List

United States of America



# Sensor Insoles

**599 \$**

*per pair*


|                  |  |
|------------------|--|
| Sensors          | 16 capacitive pressure sensors<br>6-axis inertial measurement unit |
| Sample rate      | 100/200 Hz   |
| Data acquisition | Live stream in realtime<br>to mobile device                        |
| Connectivity     | Bluetooth Low Energy 5.0   |
| Power supply     | Li-Ion PD2032 coin cell battery,<br>3.7 V, USB-C charge connector  |
| Sizes            | 9 double sizes S1 - S9<br>(EU 32 - 49, US men 1 - 13.5)            |

# Tests

*from* **0.99 \$**

*per test*

|                |  |
|----------------|--|
| Service scope  | Automated computation of test<br>results for reporting                   |
| Free services  | Test execution<br>Sensor data cloud storage<br>Biofeedback training      |
| Quota          | 10 free tests on first purchase,<br>no minimum usage,<br>no subscription |
| Tests overview | <a href="https://moticon.com/rego/tests">moticon.com/rego/tests</a>      |



All prices are net prices. Local tax applies.  
Billing period for all software as a service and cloud software services is monthly.  
Our general terms and conditions of sale apply • [moticon.com/terms-and-conditions](https://moticon.com/terms-and-conditions)  
Our terms of use and license apply • [moticon.com/terms-of-use-license](https://moticon.com/terms-of-use-license)

# Tests

---

## 1 FUNCTIONAL MOVEMENT TESTS

|        |                          |        |         |
|--------|--------------------------|--------|---------|
| RT 018 | Balance Squat            | 1 Test | 0.99 \$ |
| RT 011 | Y-Balance                | 1 Test | 0.99 \$ |
| RT 003 | Balance Front Hop        | 1 Test | 0.99 \$ |
| RT 010 | Single Leg Jump Distance | 1 Test | 0.99 \$ |
| RT 004 | Balance Side Hop         | 1 Test | 0.99 \$ |
| RT 005 | 30 Seconds Side Hop      | 1 Test | 0.99 \$ |

---

## 2 PERFORMANCE DIAGNOSTICS TESTS

|        |                       |        |         |
|--------|-----------------------|--------|---------|
| RT 013 | Counter Movement Jump | 1 Test | 0.99 \$ |
| RT 008 | Drop Jump             | 1 Test | 0.99 \$ |
| RT 037 | Foot Tapping          | 1 Test | 0.99 \$ |

---

## 3 GAIT ANALYSIS & RUNNING TESTS

|        |                               |        |         |
|--------|-------------------------------|--------|---------|
| RT 038 | Free Distance Walk and Return | 1 Test | 1.99 \$ |
| RT 039 | 15 Seconds Treadmill Walk     | 1 Test | 1.99 \$ |
| RT 040 | 15 Seconds Treadmill Run      | 1 Test | 1.99 \$ |

---

## 4 CLINICAL TESTS

|        |                          |        |         |
|--------|--------------------------|--------|---------|
| RT 001 | 2 Minute Walk            | 1 Test | 1.99 \$ |
| RT 035 | 10 Meter Walk and Return | 1 Test | 1.99 \$ |
| RT 012 | Single Leg Stance        | 1 Test | 1.99 \$ |
| RT 036 | Timed Up and Go          | 1 Test | 1.99 \$ |

All prices are net prices. Local tax applies.  
Billing period for all software as a service and cloud software services is monthly.  
Our general terms and conditions of sale apply • [moticon.com/terms-and-conditions](https://moticon.com/terms-and-conditions)  
Our terms of use and license apply • [moticon.com/terms-of-use-license](https://moticon.com/terms-of-use-license)

# Software

---

## 1 APPS

|            |  |         |
|------------|--|---------|
| Mobile App | iOS/Android, user management, test execution, test results, settings | 0.00 \$ |
| Web App    | Browser based, user management, individual & team reports, settings  | 0.00 \$ |

---

## 2 INTERFACES

|          |   |         |            |
|----------|---|---------|------------|
| Raw Data | Export of sensor raw data   | 1 month | 149.00 \$  |
| API      | Application programming interface<br>incl. raw sensor data access | 1 month | 2995.00 \$ |

All prices are net prices. Local tax applies.  
Billing period for all software as a service and cloud software services is monthly.  
Our general terms and conditions of sale apply • [moticon.com/terms-and-conditions](https://moticon.com/terms-and-conditions)  
Our terms of use and license apply • [moticon.com/terms-of-use-license](https://moticon.com/terms-of-use-license)

# Accessories & Services

## 1 POWER MANAGEMENT

|                   |  |      |          |
|-------------------|--|------|----------|
| Power Supply      | Power supply 240/110 V DC / 5 V AC, USB 2.0 type A, for 4 dual charge cables | 1 pc | 99.00 \$ |
| Dual Charge Cable | USB 2.0 type A to dual type C charge cable for 1 pair of sensor insoles      | 1 pc | 9.90 \$  |

## 2 BAGS & MARKERS

|                     |   |      |           |
|---------------------|---|------|-----------|
| Product Bag         | Cordura zipper bag for 6 pairs of sensor insoles & accessories  | 1 pc | 149.00 \$ |
| Sensorsohlensack    | Mesh sack with drawcord lock for 1 pair of sensor insoles   | 1 pc | 9.90 \$   |
| Floor Marker Disc   | Flexible floor marker with rubberized back, indicates clockwise & counter clockwise, diameter 18 cm / 7.1 in  | 1 pc | 19.90 \$  |
| Floor Marker Square | Flexible floor marker with rubberized back, indicates hop distance target, square 40 x 40 cm / 15.7 x 15.7 in | 1 pc | 29.90 \$  |
| Tape Measure        | For setting out gait corridors  | 1 pc | 9.90 \$   |

## 3 CONSUMABLES

|                     |  |         |          |
|---------------------|--|---------|----------|
| Test Socks          | Disposable socks, antibacterial, one size fits all | 100 pcs | 99.00 \$ |
| Desinfectant Spray  | Sterillium disinfectant spray, refillable          | 100 ml  | 9.90 \$  |
| Desinfectant Refill | Sterillium disinfectant, refill pack               | 1 L     | 9.90 \$  |

## 4 SERVICES

|                  |  |              |           |
|------------------|--|--------------|-----------|
| Initial Setup    | Setup customer account, account activation           | per customer | 99.00 \$  |
| Product Training | Comprehensive online training, questions and answers | 2 hr         | 299.00 \$ |

All prices are net prices. Local tax applies.  
Billing period for all software as a service and cloud software services is monthly.  
Our general terms and conditions of sale apply • [moticon.com/terms-and-conditions](https://moticon.com/terms-and-conditions)  
Our terms of use and license apply • [moticon.com/terms-of-use-license](https://moticon.com/terms-of-use-license)



Price List United States of America  
Original English Version  
March 2024

Moticon ReGo AG  
Machtlfinger Str. 21  
81379 Munich  
Germany

+49 89 200 030 160  
sales@moticon.com  
moticon.com/de/ReGo

REV 01.06