



**ReGo**  
byMOTICON

## **Pricelist**

USD • United States & Canada

## Sensor Insoles & Tests

---

### 1 **SENSOR INSOLES**

Sensors	16 capacitive pressure sensors 6-axis inertial measurement unit	per pair	\$ 949.00
Sample rate	100/200 Hz		
Data acquisition	Live stream in realtime to mobile device (iOS, Android)		
Connectivity	Bluetooth Low Energy 5.0		
Power supply	Li-Ionen PD2032 coin cell battery, 3.7 V, USB-C charge connector		
Sizes	9 double sizes S1 - S9 (EU 32 - 49, US men 1 - 13.5)		

---

### 2 **TESTS**

Service scope	Automated computation of test Results for report generation	from	\$ 0.99
Free services	Test execution Sensor data archiving Biofeedback training		
Quota	10 free tests on first purchase, no minimum usage, no subscription		
Tests overview	<a href="https://moticon.com/rego/tests">moticon.com/rego/tests</a> Test portfolio is constantly extended		

# Tests

---

## 1 **GAIT ANALYSIS & RUNNING**

RT 038	Free Distance Walk and Return	1 Test	\$ 2.99
RT 039	15 Seconds Treadmill Walk	1 Test	\$ 2.99
RT 040	15 Seconds Treadmill Run	1 Test	\$ 2.99

---

## 2 **HOP & JUMP TESTS**

RT 003	Balance Front Hop	1 Test	\$ 0.99
RT 010	Single Leg Jump Distance	1 Test	\$ 0.99
RT 004	Balance Side Hop	1 Test	\$ 0.99
RT 005	30 Seconds Side Hop	1 Test	\$ 0.99
RT 013	Counter Movement Jump	1 Test	\$ 0.99
RT 008	Drop Jump	1 Test	\$ 0.99

---

## 3 **BALANCE TESTS**

RT 012	Single Leg Stance	1 Test	\$ 0.99
RT 018	Balance Squat	1 Test	\$ 0.99
RT 011	Y-Balance	1 Test	\$ 0.99

---

## 4 **SPEED TESTS**

RT 037	Foot Tapping	1 Test	\$ 0.99
--------	--------------	--------	---------

---

## 5 **MOBILITY TESTS**

RT 001	2 Minute Walk	1 Test	\$ 2.99
RT 018	10 Meter Walk and Return	1 Test	\$ 2.99
RT 011	Timed Up and Go	1 Test	\$ 2.99

ReGo Tests are billed separately according to usage.

# Software

---

## 1 APPS

Mobile App	iOS/Android, user management, test execution, test results, settings		\$ 0.00
Web App	Browser based, user management, individual & team reports, settings		\$ 0.00

---

## 2 INTERFACES

Raw data	Export of sensor raw data	1 month	\$ 249.00
API	Application programming interface incl. raw sensor data access	1 month	\$ 1,999.00

## Accessories & Services

---

### 1 POWER MANAGEMENT

Power Supply	Power supply 110 V DC / 5 V AC, USB 2.0 Type A, for 4 dual charge cables	1 pc	\$ 99.00
Dual Charge Cable	USB 2.0 Type A dual charge cable for 1 pair of sensor insoles	1 pc	\$ 9.90

---

### 2 BAGS & MARKERS

Product Bag	Cordura zipper bag for 6 pairs of sensor insoles & accessories	1 pc	\$ 149.00
Sensor Insole Sack	Mesh sack with drawcord lock for 1 pair of sensor insoles	1 pc	\$ 9.90
Floor Marker Disc	Flexible floor marker with rubberized back, indicates clockwise & counter clockwise, diameter 18 cm / 7.1 in	1 pc	\$ 19.90
Floor Marker Square	Flexible floor maker with rubberized back, indicates hop distance target square 40 x 40 cm, 15.7 x 15.7 in	1 pc	\$ 19.90
Tape Measure	For setting out gait corridors	1 pc	\$ 9.90

---

### 3 CONSUMABLES

Test Socks	Disposable socks, antibacterial, one size fits all	100 pc	\$ 49.90
Disinfectant Spray	Sterillium disinfectant spray, refillable	100 ml	\$ 9.90
Disinfectant Refill	Sterillium disinfectant, refill pack	1 L	\$ 49.90

---

### 4 SERVICES

Initial Setup	Setup customer account, account activation	per customer	\$ 0.99
Product Training	Comprehensive online training, questions and answers	2 h	\$ 299.00

# Bundles

---

## 1 COMPLETE BUNDLE

Delivery scope	Sensor Insoles	7 pairs	\$ 7,598.20
	Power Supply	1 pc	
	Dual Charge Cable	4 pcs	
	Product Bag	1 pc	
	Floor Marker Disc	1 pc	
	Floor Marker Square	1 pc	
	Test Socks	100 pcs	
	Disinfectant Spray	100 ml	
	Initial Setup	per customer	
	Product Training	2 h	

---

## 2 STANDARD BUNDLE

Delivery scope	Sensor Insoles	5 pairs	\$ 5,620.60
	Power Supply	1 pc	
	Dual Charge Cable	2 pcs	
	Product Bag	1 pc	
	Test Socks	100 pcs	
	Disinfectant Spray	100 ml	
	Initial Setup	per customer	
	Product Training	2 h	

---

## 3 STARTER BUNDLE

Delivery scope	Sensor Insoles	2 pairs	\$ 2,733.70
	Power Supply	1 pc	
	Dual Charge Cable	1 pc	
	Product Bag	1 pc	
	Test Socks	100 pcs	
	Disinfectant Spray	100 ml	
	Initial Setup	per customer	
	Product Training	2 h	

ReGo Tests are billed separately according to usage.

All prices are net prices. Local tax applies  
Billing period is monthly for all software-as-a-service and cloud software services  
Our general terms and conditions of sale apply • [moticon.com/terms-and-conditions](https://moticon.com/terms-and-conditions)  
Our terms of usage and license apply • [moticon.com/terms-of-use-license](https://moticon.com/terms-of-use-license)



Original English version  
December 2024

Copyright  
Moticon ReGo AG  
Machtlfinger Str. 21  
81379 Munich  
Germany

+49 89 200 030 160  
[sales@moticon.com](mailto:sales@moticon.com)

REV 02.00